

# KOVACS' REAGENT

## **DISCUSSION:**

Kovacs reagent is the color indicator in the indole test. Key Kovacs Reagent is suitable for use in any classic indole procedure. For directions consult the product insert for the indole test in use.

## **STORAGE:**

Store at 2-8°C protected from light.

## **MSDS DATA:**

Kovacs reagent contains p-dimethyl- amino-benzaldehyde in Hydrochloric acid and N-Butanol. Kovacs reagent is harmful if inhaled, ingested, or absorbed through the skin. It may cause burns or severe irritation. In case of contact, immediately remove contaminated clothing and flood the skin and/or eyes with copious amounts of water. Wash contaminated garments before reusing. To the best of our knowledge, the chemical, physical and toxicological properties have not been thoroughly investigated.



KEY SCIENTIFIC PRODUCTS, INC  
1113 EAST REYNOLDS STREET  
STAMFORD, TEXAS 79553  
VOICE 800-843-1539  
FAX 888-440-4208  
WWW.KEYSCIENTIFIC.COM